

VibroTrac Table™

INTERMITTENT TRACTION



The operating principle of the VibroTrac Table is to employ vibration to evoke specific desired systemic responses (neurological and circulatory) and to utilize quadrant independent intermittent traction to evoke predictable physiological responses.

QUADRANT SPECIFIC THERAPY

- Cervical, multi-level thoracic and lumbar traction
- High to Low Intermittent traction table
- Traction force independently adjustable:
 - Right to Left
 - Head to Foot
- Adjustable table height
- Easy to use central control
- Multiple traction force options
- Specific protocols for acute/chronic, herniations, sciatica and radiculopathy
- Adjustable amplitude and frequency
- Specially designed Pelvic Belt that allows for anterior, medial, and posterior force



Patent 8257285 B2

The VibroTrac Table allows for longitudinal traction on the spine, with an emphasis on “unwinding” the spine using opposite lateral flexion to achieve de-rotation of vertebra.

DISC HYDRATION

When a disc bulges on one side there is a wedging on the opposite side. The VibroTrac Table protocol maximizes traction on the wedged anterior side. This maximizes the billows effect, increases circulation, reduces swelling and promotes resorption of nuclear material there by reducing impingement.



Knowles, J Pre 1/16/08

Knowles, J Post 10/22/08

Knowles, J Post 1/22/15

“After three weeks of numbness in my low back and foot from a herniated disc, I was scheduled to have low back surgery. An MRI showed that I had two herniated discs. After one treatment utilizing the PneuBack Program, I cancelled my scheduled surgery. The numbness in my right leg is now 95% gone after 8 visits and the disc has gotten better. I feel NO PAIN!”

Chet Howell

“When I started the PneuBack Program, I was having back pain that prevented me from going about my daily activities. I was looking at the possibility of surgery. The PneuBack Program gave me another option and took the pain away immediately. I now have a long-term plan and feel that we have treated the cause and not just the symptoms.”

Genuine Pneumex Equipment

www.pneumex.com

pneumex@pneumex.com

800.447.5792

Specifications

VibroTrac™

INTERMITTENT TRACTION

ACCESSORIES INCLUDED:

Compressor


Specially made Oil-less to run quieter In your clinic.



Table Vests

X-Heart


One Small
One Medium
One Large



Ball 55cm



VibroTrac Wedges



4" & 6"

Ankle Tie



Mobilization Belt



72"

Mobilization Pad



2"

Pelvic Belt

One Small
One Medium
One Large



Pneumex IRB Study

IBI (PneuBack Program Provider) / FOTO

Published in peer review journal on line 8/26/2016

www.journalchiromed.com

Mean number of visitations

FOTO: 13.43

IBI: 13.60

Pain scale results (0-10)

FOTO: First Visit: 7.978

Second Visit: 5.532

IBI: First Visit: 7.474

Second Visit: 2.474

Oswestry

FOTO: First Visit: 41.83

Second Visit: 32.72

IBI: First Visit: 46.33

Second Visit: 25.513

Mean change score for Pain

FOTO: 2.4

IBI: 5.0

Mean change score for Oswestry

FOTO: 9.11

IBI: 20.82

SPECIFICATIONS:

Height	Dimensions
18 to 44"	30x 78"
Capacity	Speed
400 lbs	10 to 60 Hz

DISTRIBUTED BY:

Pneumex, Inc.

2605 North Boyer Ave
Sandpoint, Idaho 83864
800.447.5792



Patented Innovation
For Rehabilitation

Genuine Pneumex Equipment

www.pneumex.com

pneumex@pneumex.com

800.447.5792

Rev 2-16-2017