

PNEUVEST™

3 IN 1 UNWEIGHTING HARNESS



1. Using both the upper and the lower belt the weight removed is dispersed evenly between the waist and the buttock.

2. Using only the upper belt a patient can be lowered onto a commode or positioned to change pants etc. without interference from the lower strap.

3. If the patient is unable to take pressure, the upper belt can be undone so that most of the weight is dispersed onto the lower belt.

MULTI-PURPOSE HARNESS

- Worry-free unweighted activities.
- Exercise in an upright, controlled and safe position.
- Comfortable, normal rib and diaphragm interaction.
- Doesn't ride up.
- Can be 100% Upper or 100% Lower.
- Unisex design fits both male and female.
- Easy wheelchair use.



Upper and Lower Belt

1

Patent Pending

Upper Belt

2



Lower Belt

3



Description	Waist Size	Color
Small	22" to 32"	Red
Medium	32" to 44"	Blue
Large	44" to 55"	Yellow

"Every time my back went out I would spend days in bed. I tried everything short of surgery. The clinic in Sandpoint introduced me to the PneuBack Program. I am doing things I thought I would never do again. I ran with my dog for the first time in 15 years."

Wes

"We love the PneuBack Program! We see a lot of back injury patients and most of them have already been to other clinics over their years of back pain. Once started in the program, patients often wish they could have come years ago."

Don Gale, PT

Genuine Pneumex Equipment

www.pneumex.com

info@pneumex.com

800.447.5792